



# Labcor

*pediatrics*

[www.labcor.com.br](http://www.labcor.com.br)

 **Labcor**  
*Commitment to life*

# VENA<sup>®</sup>

*Graft Cardiovascular Biology*

Biological cardiovascular graft from bovine jugular vein, treated with glutaraldehyde. It functions as a valved vascular substitute for veins stricken by diverse pathologies, particularly atresia or stenosis. A one piece unit in which only the collaterals are linked.

## Available Size

Ø 11\*, 13, 15, 17 e 19 mm x 100 mm  
(\*Available thru special order)

Vena<sup>®</sup> is indicated for extracardiac Fontan surgery.

# VENA<sup>®</sup>

*Enxerto Cardiovascular Biológico*

Enxerto cardiovascular biológico obtido a partir da veia jugular bovina, tratado em glutaraldeído. Funciona como substituto cardiovascular valvado para vasos acometidos por patologias diversas, especialmente atresia ou estenose. Apresenta-se íntegra e apenas as colaterais são ligadas.

## Available Size

Ø 11\* , 13, 15, 17 e 19 mm x 100 mm  
(\*Consulte o representante comercial)

Vena<sup>®</sup> é indicado para cirurgias de Fontan extracardíaco.



Av. Silva Lobo, 1670, Blocos 1 e 2, Nova Granada  
Belo Horizonte, MG - Brasil - CEP 30.431-262  
Phone: + 55 31 3371 1371 \*contato@labcor.com.br

[www.labcor.com.br](http://www.labcor.com.br)

P/N: 99B0141A