

TO WHOM IT MAY CONCERN

Labcor informs that the products listed below are MR safe. Specific materials used to make those products are essentially nonmagnetic, nonmetallic and non-conducting.

Name	Type	Sizes
MAGIS	Annuloplasty Ring	ALL Sizes
T07 EAIS	Biological Conduit	ALL Sizes
VELUM L-HYDRO	Biological Conduit	ALL Sizes
VASTI L-HYDRO	Biological Conduit	ALL Sizes
VIVENDI L-HYDRO	Biological Patch	ALL Sizes
BABYGRAFT/ COROGRAFT	Biological Graft	ALL Sizes
JP / LAB	Biological Monocuspid Patch	ALL Sizes
VENA	Biological Conduit	ALL Sizes
PERICADIUM PATCH	Biological Patch	ALL Sizes