

## TO WHOM IT MAY CONCERN

MRI testing has demonstrated that the product Star – Valvular annuloplasty ring – Model Tricuspid is MR Conditional.

A patient with the device Star – Valvular annuloplasty ring – Model Tricuspid can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 10,000-Gauss/cm (extrapolated) or less
- Maximum MR system reported, whole body averaged specific absorption rate of 4-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the First Level Controlled Operating Mode of operation for the MR system.

<b>Name</b>	<b>Type</b>	<b>Sizes</b>	<b>MR</b>
Star	Annuloplasty Ring	All sizes	MR Conditional